

Aquaporin 4 Rabbit pAb


ABZ1376

Product Overview


Name	Aquaporin 4 Rabbit pAb
Accession(Primary)	P55087
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (P55087)

Gene Name	AQP4
About Gene Symbol	



Western blot analysis of AQP4 in mouse heart lysates using AQP4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using Aquaporin 4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.