

FGF8 Rabbit pAb

ABZ1377

Product Overview

Name	FGF8 Rabbit pAb
Accession(Primary)	P55075
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (P55075)

Gene Name FGF8**About Gene Symbol**

Western blot analysis of FGF8 in K562 lysates
using FGF8 antibody.

