

LIM Kinase 1 Rabbit pAb


ABZ1379

Product Overview

Name	LIM Kinase 1 Rabbit pAb
Accession(Primary)	P53667
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P53667)

Gene Name	LIMK1
About Gene Symbol	



Western blot analysis of LIM Kinase 1 in mouse brain lysates using LIM Kinase 1 antibody.