

JAK3 Rabbit pAb


ABZ1380

Product Overview


Name	JAK3 Rabbit pAb
Accession(Primary)	P52333
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 125 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P52333)


Gene Name	JAK3
About Gene Symbol	




Immunofluorescence analysis of JAK3 in Human uterus tissue using JAK3 antibody(red),and DAPI (blue).




Immunofluorescence analysis of JAK3 in rat kidney using JAK3 antibody(red),and DAPI (blue).




Immunohistochemistry analysis of paraffin-embedded Human colon tissue using JAK3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded Human liver tissue using JAK3 antibody.High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded Human stomach tissue using JAK3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunohistochemistry analysis of paraffin-embedded Human stomach cancer tissue using JAK3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.