

CCR5 Rabbit pAb

ABZ1381


Product Overview

Name	CCR5 Rabbit pAb
Accession(Primary)	P51681
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (P51681)

Gene Name CCR5

About Gene Symbol



Western blot analysis of CCR5 in various lysates using CCR5 antibody.

