

VASP Rabbit pAb


ABZ1382

Product Overview


Name	VASP Rabbit pAb
Accession(Primary)	P50552
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 46,50 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P50552)

Gene Name	VASP
About Gene Symbol	VASP (Vasodilator Stimulated Phosphoprotein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of VASP in various lysates using VASP antibody.



Western blot analysis of VASP in NIH/3T3 lysates treated with forskolin using VASP antibody. The lane on the right is blocked with the synthesized peptide.