

SDF1 Rabbit pAb


ABZ1386

Product Overview

Name	SDF1 Rabbit pAb
Accession(Primary)	P48061
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (P48061)

Gene Name	CXCL12
About Gene Symbol	



Western blot analysis of SDF1 in various lysates using SDF1 antibody.