



Melatonin Receptor 1A Rabbit pAb

ABZ1387

Product Overview

Name	Melatonin Receptor 1A Rabbit pAb
Accession(Primary)	P48039
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300; IF-1:50-1:200

Target information (P48039)

Gene Name MTNR1A**About Gene Symbol**

Immunofluorescence analysis of Melatonin Receptor 1A in HepG2 cells using Melatonin Receptor 1A antibody. Sample with blocking peptide on the right.