

Fas Ligand Rabbit pAb


ABZ1388

Product Overview

Name	Fas Ligand Rabbit pAb
Accession(Primary)	P48023
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Pig
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:40000

Target information (P48023)

Gene Name	FASLG
About Gene Symbol	



Western blot analysis of Fas Ligand in various lysates using Fas Ligand antibody.