

ASH2L Rabbit mAb

ABZ1389

Product Overview



Name	ASH2L Rabbit mAb
Accession(Primary)	Q9UBL3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20-1:50

Target information (Q9UBL3)

Gene Name	ASH2L
About Gene Symbol	



Immunocytochemistry analysis of ASH2L (green) in K562 using ASH2L antibody, and DAPI(blue).



Western blot analysis of ASH2L in K562, C6 lysates using ASH2L antibody.