

Syk Rabbit pAb

ABZ1390


Product Overview

Name	Syk Rabbit pAb
Accession(Primary)	P43405
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (P43405)

Gene Name SYK**About Gene Symbol**

Western blot analysis of Tyrosine Protein Kinase SYK in Raji lysates using Tyrosine Protein Kinase SYK antibody



Immunohistochemistry analysis of paraffin-embedded Human lymphoma using Tyrosine Protein Kinase SYK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.