


CXCL5 Rabbit pAb

ABZ1391


Product Overview

Name	CXCL5 Rabbit pAb
Accession(Primary)	P42830
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200

Target information (P42830)

Gene Name CXCL5**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human pancreascancer using ENA 78 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using ENA 78 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.