

Bcl 6 Rabbit pAb

ABZ1392

Product Overview

Name	Bcl 6 Rabbit pAb
Accession(Primary)	P41182
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 79 kDa; Observed MW: 79 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P41182)

Gene Name	BCL6
About Gene Symbol	



Western blot analysis of Bcl6 in 293 lysates using Bcl6 antibody.