

STAT3 Rabbit pAb


ABZ1393

Product Overview


Name	STAT3 Rabbit pAb
Accession(Primary)	P40763
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (P40763)


Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of STAT3 in various lysates using STAT3 antibody.



Western blot analysis of STAT3 in NIH3T3 lysates using Stat3 antibody.



Western blot analysis of STAT3 in Hepg2 lysates using STAT3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using STAT3 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was preabsorbed by immunogen peptide.