

CD130 Rabbit pAb


ABZ1394

Product Overview


Name	CD130 Rabbit pAb
Accession(Primary)	P40189
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 104 kDa; Observed MW: 130-140kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P40189)

Gene Name	IL6ST
About Gene Symbol	



Western blot analysis of CD130 in 293T, COLO lysates using CD130 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain using CD130 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.