


p21 Rabbit pAb





ABZ1395


Product Overview

Name	p21 Rabbit pAb
Accession(Primary)	P38936
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:50-1:200; IHC-1:50-1:100

Target information (P38936)

Gene Name CDKN1A**About Gene Symbol**





Western blot analysis of p21 in various lysates
using p21 antibody