

PPAR gamma Rabbit pAb


ABZ1397

Product Overview


Name	PPAR gamma Rabbit pAb
Accession(Primary)	P37231
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 53,57 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P37231)

Gene Name	PPARG
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse Heart Tissue using PPAR gamma antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PPAR gamma in HeLa, 3T3, PC-12 lysates using PPAR gamma antibody