

TGF beta Receptor II Rabbit pAb


ABZ1399

Product Overview

Name	TGF beta Receptor II Rabbit pAb
Accession(Primary)	P37173
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P37173)

Gene Name	TGFBR2
About Gene Symbol	



Western blot analysis of TGF beta Receptor II in Jurkat lysates using TGF beta Receptor II antibody.