

SREBP1 Rabbit pAb


ABZ1400

Product Overview

Name	SREBP1 Rabbit pAb
Accession(Primary)	P36956
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 122 kDa; Observed MW: 122 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P36956)

Gene Name	SREBF1
About Gene Symbol	



Western blot analysis of SREBP1 in HeLa, HepG2 lysates using SREBP1 antibody.