

CD135 Rabbit pAb


ABZ1401

Product Overview

Name	CD135 Rabbit pAb
Accession(Primary)	P36888
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 113 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (P36888)

Gene Name	FLT3
About Gene Symbol	



Western blot analysis of CD135 in 4T1 lysates using CD135 antibody.