

TIMP3 Rabbit pAb

ABZ1403


Product Overview

Name	TIMP3 Rabbit pAb
Accession(Primary)	P35625
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (P35625)

Gene Name TIMP3**About Gene Symbol**

Western blot analysis of Metalloproteinase
Inhibitor 3 in Huvec lysates using Metalloproteinase
Inhibitor 3 antibody



Immunohistochemistry analysis of paraffin-embedded Human cervical cancer using Metalloproteinase Inhibitor 3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.