

Thrombospondin 2 Rabbit pAb


ABZ1404

Product Overview


Name	Thrombospondin 2 Rabbit pAb
Accession(Primary)	P35442
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 130 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:50-1:200; ELISA-1:10000; IHC-1:50-1:100

Target information (P35442)

Gene Name	THBS2
About Gene Symbol	



Western blot analysis of Thrombospondin 2 in mouse heart lysates using Thrombospondin 2 antibody.



Immunohistochemistry analysis of paraffin-embedded rat testis using Thrombospondin 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.