

VE Cadherin Rabbit pAb


ABZ1405

Product Overview


Name	VE Cadherin Rabbit pAb
Accession(Primary)	P33151
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 88 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (P33151)

Gene Name	CDH5
About Gene Symbol	



Western blot analysis of Cadherin 5 in Jurkat lysates using Cadherin 5 antibody.



Western blot analysis of Cadherin 5 in A549 lysates using Cadherin 5 antibody.