

14-3-3 Rabbit pAb


ABZ1406

Product Overview


Name	14-3-3 Rabbit pAb
Accession(Primary)	P31946
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27-29 kDa; Observed MW: 27-29 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P31946)

Gene Name	YWHAB/YWHAЕ/YWHAH
About Gene Symbol	YWHAB is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of 14-3-3 (pan) in rat Hippocampus, rat Cerebral Cortex and rat Brain lysates using 14-3-3 (pan) antibody.



Western blot analysis of 14-3-3 in C6, 3T3 and HeLa lysates using 14-3-3 (pan) antibody.