

ASH2L Rabbit mAb


ABZ1407

Product Overview


Name	ASH2L Rabbit mAb
Accession(Primary)	Q9UBL3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:200-1:1000

Target information (Q9UBL3)


Gene Name	ASH2L
About Gene Symbol	



Western blot analysis of ASH2L in C6, 3T3, HeLa lysates using ASH2L antibody.



Immunocytochemistry analysis of ASH2L(green) in HeLa using ASH2L antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ASH2L antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.