

## AKT1/2/3 Rabbit pAb


ABZ1408

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | AKT1/2/3 Rabbit pAb                        |
| <b>Accession(Primary)</b>   | P31749                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 56 kDa; Observed MW: 60 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (P31749)

|                          |                                                                                                                                                                                                                                                                        |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | AKT1/AKT2/AKT3                                                                                                                                                                                                                                                         |
| <b>About Gene Symbol</b> | AKT1 (AKT Serine/Threonine Kinase 1) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of AKT in 3T3, C6 and HeLa lysates using Akt antibody.