

CDC25A Rabbit pAb


ABZ1410

Product Overview


Name	CDC25A Rabbit pAb
Accession(Primary)	P30304
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P30304)

Gene Name	CDC25A
About Gene Symbol	



Western blot analysis of CDC25A in 231, A549 lysates using CDC25A antibody



Immunohistochemistry analysis of paraffin-embedded Human thyroid cancer using CDC25A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.