

eNOS Rabbit pAb

ABZ1411


Product Overview

Name	eNOS Rabbit pAb
Accession(Primary)	P29474
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 133 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:50-1:300; IHC-1:50-1:100


Target information (P29474)

Gene Name NOS3


About Gene Symbol




Western blot analysis of eNOS in HepG2 lysates using eNOS antibody.



Western blot analysis of eNOS in HepG2 lysates using NOS3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using eNOS antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of eNOS in Jurkat lysates treated with Insulin using eNOS antibody. The lane on the right is blocked with the synthesized peptide.