

Caspase 1 Rabbit pAb

ABZ1412

Product Overview


Name	Caspase 1 Rabbit pAb
Accession(Primary)	P29466
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300; IF-1:50-1:200; ELISA-1:10000-1:20000

Target information (P29466)

Gene Name	CASP1
About Gene Symbol	



Western blot analysis of Caspase 1 in various lysates using Caspase 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human stomach using Caspase 1 antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval.