

NMBR Rabbit pAb


ABZ1414

Product Overview

Name	NMBR Rabbit pAb
Accession(Primary)	P28336
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-P-1:50-1:200; IF-1:200-1:1000

Target information (P28336)

Gene Name	NMBR
About Gene Symbol	



Western blot analysis of Neuromedin B Receptor
in Human Brain, rat Brain lysates using
NMBR(Bombesin Receptor)antibody.