

Calnexin Rabbit pAb


ABZ1415

Product Overview

Name	Calnexin Rabbit pAb
Accession(Primary)	P27824
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P27824)

Gene Name	CANX
About Gene Symbol	CANX (Calnexin) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Calnexin in CHO-K1, Jurkat, 3T3, C6 and HeLa lysates using Calnexin antibody.