

IGFBP3 Rabbit pAb


ABZ1416

Product Overview


Name	IGFBP3 Rabbit pAb
Accession(Primary)	P17936
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:500; IF/ICC1:100-1:500; ELISA-1:5000-1:20000

Target information (P17936)


Gene Name	IGFBP3
About Gene Symbol	




Western blot analysis of IGFBP3 in HeLa, HepG2 lysates using IGFBP3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using IGFBP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human kidney using IGFBP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of IGFBP3 in mosue brian
lysates using IGFBP3 antibody.