

Desmin Rabbit pAb


ABZ1417

Product Overview


Name	Desmin Rabbit pAb
Accession(Primary)	P17661
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P17661)


Gene Name	DES
About Gene Symbol	DES (Desmin) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cardiovascular disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Desmin in various lysates using Desmin antibody.




Western blot analysis of Desmin in VEC lysates using Desmin antibody.




Western blot analysis of Desmin in pancreatic stellate cell lysates using Desmin antibody.

Union Hospital of Wuhan



Immunohistochemistry analysis of paraffin-embedded Human skeletal muscle tissue using Desmin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of Desmin in K562 lysates using Desmin antibody. The lane on the right is blocked with the synthesized peptide.