

## PKA alpha/beta/gamma Rabbit pAb


ABZ1418

### Product Overview

|                             |                                                                     |
|-----------------------------|---------------------------------------------------------------------|
| <b>Name</b>                 | PKA alpha/beta/gamma Rabbit pAb                                     |
| <b>Accession(Primary)</b>   | P17612                                                              |
| <b>Host Species</b>         | Rabbit                                                              |
| <b>Clonality</b>            | Polyclonal                                                          |
| <b>Isotype</b>              | IgG                                                                 |
| <b>Species Reactivity</b>   | Human, Mouse, Rat, Pig                                              |
| <b>Calculated MW</b>        | Calculated MW: 40 kDa; Observed MW: 40 kDa                          |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:200-1:1000; ELISA-1:20000 |

### Target information (P17612)

|                          |                                                                                                                                                                                                                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | PRKACA/PRKACB                                                                                                                                                                                                                                                                                            |
| <b>About Gene Symbol</b> | PRKACA (Protein Kinase) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of cAMP Protein Kinase Catalytic Subunit in various lysates using cAMP Protein Kinase Catalytic Subunit antibody.

Western blot analysis of cAMP Protein Kinase Catalytic Subunit in COLO205 lysates using PKA $\alpha$ / $\beta$ / $\gamma$  cat antibody.