

Connexin 43 / GJA1 Rabbit pAb


ABZ1420

Product Overview

Name	Connexin 43 / GJA1 Rabbit pAb
Accession(Primary)	P17302
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P17302)

Gene Name	GJA1
About Gene Symbol	



Western blot analysis of Connexin 43 in various
lysates using Connexin 43 antibody.