

Histone H1 Rabbit pAb


ABZ1421

Product Overview


Name	Histone H1 Rabbit pAb
Accession(Primary)	P16401
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 22 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P16401)

Gene Name	HIST1H1B
About Gene Symbol	



Western blot analysis of Histone H1 in various lysates using Histone H1 antibody.



Western blot analysis of Histone H1 in 3T3 lysates using Histone H1 antibody.