

Wnt1 Rabbit pAb


ABZ1423

Product Overview


Name	Wnt1 Rabbit pAb
Accession(Primary)	P04628
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 41 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P04628)

Gene Name	WNT1
About Gene Symbol	



Western blot analysis of Wnt1 in various lysates using Wnt1 antibody.



Western blot analysis of Wnt1 in HT-29 lysates using Wnt1 antibody.