

SOD2 Rabbit pAb


ABZ1424

Product Overview


Name	SOD2 Rabbit pAb
Accession(Primary)	P04179
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P04179)


Gene Name	SOD2
About Gene Symbol	SOD2 (Superoxide Dismutase 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of SOD2 in 293, COS7, C6, 3T3 and HeLa lysates using SOD2 antibody.



Western blot analysis of SOD2 in HepG2 and Raw264.7 lysates using SOD2 antibody.



Western blot analysis of SOD2 in K562, 293T, 3T3 and Jurkat lysates using SOD2 antibody.