

## FN Rabbit pAb


ABZ1427

### Product Overview


|                             |                                                                     |
|-----------------------------|---------------------------------------------------------------------|
| <b>Name</b>                 | FN Rabbit pAb                                                       |
| <b>Accession(Primary)</b>   | P02751                                                              |
| <b>Host Species</b>         | Rabbit                                                              |
| <b>Clonality</b>            | Polyclonal                                                          |
| <b>Isotype</b>              | IgG                                                                 |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                   |
| <b>Calculated MW</b>        | Calculated MW: 272 kDa; Observed MW: 260 kDa                        |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC:-1:100-1:300;<br>IF-1:100-1:300; ELISA-1:20000 |

### Target information (P02751)


|                          |                                                                                                                                                                                                                                                                                                 |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | FN1                                                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | FN1 (Fibronectin 1) is a protein kinase and signaling regulator involved in cell migration and invasion. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |




Western blot analysis of Fibronectin in various lysates using Fibronectin antibody.




Western blot analysis of Fibronectin in HeLa lysates using FN1 antibody.



Western blot analysis of Fibronectin in 3T3 KB K562 293T lysates using FN1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human thyroid gland tissue using Fibronectin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of Fibronectins in HeLa lysates using Fibronectin 1 antibody. The lane on the right is blocked with the synthesized peptide.