

Collagen IV alpha 1 Rabbit pAb

ABZ1429


Product Overview

Name	Collagen IV alpha 1 Rabbit pAb
Accession(Primary)	P02462
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 161 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (P02462)

Gene Name COL4A1

About Gene Symbol



Western blot analysis of Collagen IV alpha 1 in Hela MCF-7 U2OS 293T KB 293T VEC lysates using Collagen IV alpha 1 antibody.

