

WDR77 (8A10) Mouse mAb


ABZ1430

Product Overview


Name	WDR77 (8A10) Mouse mAb
Accession(Primary)	Q9BQA1
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (Q9BQA1)

Gene Name	WDR77
About Gene Symbol	WDR77 (WD Repeat Domain 77) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of WDR77 in C6, 3T3 and K562 lysates using WDR77 antibody.



Immunocytochemistry analysis of WDR77 (8A10) in HeLa using WDR77 antibody.