

IL-1 beta Rabbit pAb


ABZ1431

Product Overview


Name	IL-1 beta Rabbit pAb
Accession(Primary)	P01584
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 17,31 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P01584)


Gene Name	IL1B
About Gene Symbol	




Western blot analysis of IL-1 beta in various cell lines lysates using IL-1B antibody.




Western blot analysis of IL-1 beta in THP-1 lysates using IL-1B antibody.



Immunohistochemistry analysis of paraffin-embedded rat lung using IL-1B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of IL-1 beta in THP-1 cells using IL-1B antibody, and DAPI (blue).



Immunofluorescence analysis of IL-1 beta in mouse lung using IL-1B antibody, and DAPI (blue).