

TNF alpha Rabbit pAb


ABZ1434

Product Overview

Name	TNF alpha Rabbit pAb
Accession(Primary)	P01375
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 15-30 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P01375)

Gene Name	TNF
About Gene Symbol	



Western blot analysis of TNF alpha in Jurkat
lysates using TNF alpha antibody.