

NGF Rabbit pAb


ABZ1435

Product Overview

Name	NGF Rabbit pAb
Accession(Primary)	P01138
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P01138)

Gene Name	NGF
About Gene Symbol	



Western blot analysis of NGF in various lysates using NGF antibody.



Western blot analysis of NGF in Jurkat lysates using NGF antibody.