

TGF beta 1 Rabbit pAb


ABZ1436

Product Overview


Name	TGF beta 1 Rabbit pAb
Accession(Primary)	P01137
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P01137)


Gene Name	TGFB1
About Gene Symbol	TGFB1 (Transforming Growth Factor Beta 1) is a formin actin nucleator involved in TGF- β signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of TGF beta 1 in Jurkat and HeLa lysates using TGF beta 1 antibody.




Western blot analysis of TGF beta 1 in MCF-7 lysates using TGFbeta1 antibody.



Western blot analysis of TGF beta 1 in mouse brain lysates using TGFu3b21 antibody

mouse-brain



Western blot analysis of TGF beta 1 in mouse brain lysates using TGF beta 1 antibody

mouse-brain

