

PDGF B Rabbit pAb

ABZ1437

Product Overview

Name	PDGF B Rabbit pAb
Accession(Primary)	P01127
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P01127)

Gene Name	PDGFB
About Gene Symbol	