

GTPase NRAS+HRAS+KRAS Rabbit pAb


ABZ1438

Product Overview

Name	GTPase NRAS+HRAS+KRAS Rabbit pAb
Accession(Primary)	P01111
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (P01111)

Gene Name	NRAS
About Gene Symbol	



Western blot analysis of N/H/K-Ras in T47D, 3T3, COS7 and HeLa lysates using N/H/K-Ras antibody.