

c-Myc Rabbit pAb


ABZ1439

Product Overview


Name	c-Myc Rabbit pAb
Accession(Primary)	P01106
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 50-60 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:40000

Target information (P01106)

Gene Name	MYC
About Gene Symbol	MYC (MYC Proto-Oncogene, BHLH Transcription Factor) is a proto-oncogene involved in translational regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of c-Myc in HuvEc lysates using c-Myc antibody.



Immunohistochemistry analysis of paraffin-embedded mouse liver tissue using c-Myc antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.