

c-Fos Rabbit pAb


ABZ1440

Product Overview

Name	c-Fos Rabbit pAb
Accession(Primary)	P01100
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (P01100)

Gene Name	FOS
About Gene Symbol	



Western blot analysis of c-Fos in various lysates using c-Fos antibody.