

TIMP1 Rabbit pAb

ABZ1441


Product Overview

Name	TIMP1 Rabbit pAb
Accession(Primary)	P01033
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Rabbit
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:50-1:200; IHC-1:50-1:100

Target information (P01033)

Gene Name TIMP1

About Gene Symbol



Western blot analysis of Metalloproteinase Inhibitor 1 in various lysates using Metalloproteinase Inhibitor 1 antibody.

